

OTICON | Real

Technical data sheet

miniRITE R

60 85 100 105



	Real 1	Real 2	Real 3
Speech Understanding	MoreSound Intelligence™ 2.0	Level 1	Level 2
	- Environment configuration	5 Options	5 Options
	- Virtual Outer Ear	3 Configurations	1 Configuration
	- Spatial Balancer	100%	60%
	- Neural Noise Suppression, Difficult / Easy	10 dB / 4 dB	6 dB / 2 dB
	- Sound Enhancer	3 Configurations	2 Configurations
	- Wind & Handling Stabilizer	•	•
	MoreSound Amplifier™ 2.0	•	•
	- SuddenSound Stabilizer	6 Configurations	4 Configurations
	Feedback Prevention	MoreSound Optimizer™ & Feedback shield	MoreSound Optimizer™ & Feedback shield
Sound Quality	Spatial Sound™	4 Estimators	2 Estimators
	Soft Speech Booster	•	•
	Frequency lowering	Speech Rescue™	Speech Rescue™
	Clear Dynamics	•	-
	Better-Ear Priority	•	-
	Fitting Bandwidth ¹	10 kHz	8 kHz
	Bass Boost (streaming)	•	•
	Processing Channels	64	48
	Fitting Bands	24	18
	Multiple Directionality options	•	•
Personalization & Optimizing Fitting	Adaptation Management	•	•
	Fitting Formulas	VAC+, NAL-NL1/ NAL-NL2, DSL v5	VAC+, NAL-NL1/ NAL-NL2, DSL v5
	Connecting to the world	•	•
	Hands-free communication ²	•	•
	Direct streaming ³	•	•
	ConnectClip	•	•
	EduMic	•	•
	Remote Control 3.0	•	•
	TV Adapter 3.0	•	•
	Phone Adapter 2.0	•	•
Connecting to the world	Tinnitus SoundSupport™	•	•
	CROS/BiCROS support	•	•

1) Bandwidth accessible for gain adjustments during fitting

2) Hands-free communication is available with iPhone 11 or later running iOS 15.2 or later, and iPad running iPadOS® 15.2 or later

3) From select iPhone, iPad, iPod touch, and select Android devices with the Audio Streaming for Hearing Aids (ASHA) protocol

Operating and charging conditions

Temperature: +5°C to +40°C (41°F to 104°F)

Humidity: 5% to 93% relative humidity,

non-condensing

Atmospheric pressure: 700 hPa to 1060 hPa

Storage and transportation conditions

Temperature and humidity shall not exceed the below limits for extended periods during transportation and storage.

Transport

Temperature: -20°C to +60°C (-4°F to 140°F)

Humidity: 5% to 93% relative humidity,

non-condensing

Atmospheric pressure: 700 hPa to 1060 hPa

Storage

Temperature: -20°C to +30°C (-4°F to 86°F)

Humidity: 5% to 93% relative humidity,

non-condensing

Atmospheric pressure: 700 hPa to 1060 hPa

Apple, the Apple logo, iPhone®, iPad®, and iPod touch® are trademarks of Apple Inc., registered in the U.S. and other countries.

Oticon Real™ miniRITE R offers a discreet design. It is powered by a rechargeable lithium-ion battery and features telecoil and a double push-button. Based on Bluetooth® Low Energy technology, it is a Made for iPhone® hearing aid and supports hands-free communication and direct streaming for iPhone®, iPad®, iPod touch® and select Android™ devices.

MoreSound Intelligence™ creates a more precise and natural representation of individual sounds with clearer and more distinct contrasts.

Oticon Real is built on the Polaris R™ platform, which utilizes faster detectors for powering new innovations used to optimize the audibility of the environmental sounds in the sound scene.




For information on compatibility, please visit www.oticon.com/support/compatibility

oticon
life-changing technology

		Ear Simulator	2CC Coupler
		Measured according to IEC 60118-0:1983/AMD1:1994, IEC 60118-0:2015, IEC 60118-1:1995+AMD1:1998 CSV and IEC 60318-4:2010	Measured according to ANSI S3.22-2014, IEC 60118-0:2015 and IEC 60318-5:2006
<			

1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+AMD1:1994 but without influence of feedback.

2) Expected operating time for rechargeable battery depends on use pattern, active feature set, hearing loss, sound environment, battery age and use of wireless accessories.










1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+AMD1:1994 but without influence of feedback.

2) Expected operating time for rechargeable battery depends on use pattern, active feature set, hearing loss, sound environment, battery age and use of wireless accessories.

		Ear Simulator Measured according to IEC 60118-0:1983/AMD1:1994, IEC 60118-0:2015, IEC 60118-1:1995+AMD1:1998 CSV and IEC 60318-4:2010	2CC Coupler Measured according to ANSI S3.22-2014, IEC 60118-0:2015 and IEC 60318-5:2006
85 Mold, Bass & Power dome OpenBass dome Technical information Omnidirectional mode is used unless otherwise stated.		OSPL90 dB Full-on Gain Frequency Response	
OSPL90	Peak (dB SPL)	127	117
	1600 Hz (dB SPL)	121	113
	HFA-OSPL90 (dB SPL)	122	114
Full-on Gain ¹	Peak (dB)	66	55
	1600 Hz (dB)	53	45
	HFA-FOG (dB)	56	48
Reference test gain (dB)	46	37	
Frequency range (Hz)	100-9500	100-8900	
Telecoil output	1 mA/m field (1600 Hz) (dB SPL)	84	
	10 mA/m field (1600 Hz) (dB SPL)	104	
	HFA-SPLITS L/R (dB SPL)		94/94
Total harmonic distortion (Input 70 dB SPL)	500 Hz (%)	<2	<2
	800 Hz (%)	<4	<2
	1600 Hz (%)	<5	<2
Equivalent input noise level	Omni (dB SPL)	21	18
	Dir (dB SPL)	29	28
Battery	Lithium-ion	Lithium-ion	
Expected operating time, hours ²	24		








1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+AMD1:1994 but without influence of feedback.

2) Expected operating time for rechargeable battery depends on use pattern, active feature set, hearing loss, sound environment, battery age and use of wireless accessories.

		Ear Simulator	2CC Coupler	
		Measured according to IEC 60118-0:1983/AMD1:1994, IEC 60118-0:2015, IEC 60118-1:1995+AMD1:1998 CSV and IEC 60318-4:2010	Measured according to ANSI S3.22-2014, IEC 60118-0:2015 and IEC 60318-5:2006	
 85  — Acoustic input: 60 dB SPL - - - Magnetic input: 31.6 mA/m		OSPL90  dB  dB SPL 	OSPL90  dB  Frequency Response 	
OSPL90		Peak (dB SPL)	127	117
		1600 Hz (dB SPL)	121	113
		HFA-OSPL90 (dB SPL)	122	114
Full-on Gain ¹		Peak (dB)	66	55
		1600 Hz (dB)	53	45
		HFA-FOG (dB)	56	48
Reference test gain (dB)			46	37
Frequency range (Hz)			100-7500	100-7500
Telecoil output		1 mA/m field (1600 Hz) (dB SPL)	84	
		10 mA/m field (1600 Hz) (dB SPL)	104	
		HFA-SPLITS L/R (dB SPL)		94/94
Total harmonic distortion (Input 70 dB SPL)		500 Hz (%)	<2	<2
		800 Hz (%)	<4	<2
		1600 Hz (%)	<5	<2
Equivalent input noise level		Omni (dB SPL)	22	18
		Dir (dB SPL)	29	27
Battery			Lithium-ion	Lithium-ion
Expected operating time, hours ²			24	








1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+AMD1:1994 but without influence of feedback.

2) Expected operating time for rechargeable battery depends on use pattern, active feature set, hearing loss, sound environment, battery age and use of wireless accessories.

		Ear Simulator	2CC Coupler	
		Measured according to IEC 60118-0:1983/AMD1:1994, IEC 60118-0:2015, IEC 60118-1:1995+AMD1:1998 CSV and IEC 60318-4:2010	Measured according to ANSI S3.22-2014, IEC 60118-0:2015 and IEC 60318-5:2006	
 100 Power receiver mold, Bass & Power dome Technical information Omnidirectional mode is used unless otherwise stated. Warning to the hearing aid dispenser The maximum output capability of the hearing aid may exceed 132 dB SPL (IEC 711). Special care should be exercised in selecting and fitting the hearing aid, as there may be risk of impairing the remaining hearing of the hearing aid user. Acoustic input: 60 dB SPL Magnetic input: 31.6 mA/m		OSPL90 	OSPL90  Full-on Gain  Frequency Response 	
OSPL90		Peak (dB SPL)	132	124
		1600 Hz (dB SPL)	130	122
		HFA-OSPL90 (dB SPL)	127	120
Full-on Gain ¹		Peak (dB)	66	57
		1600 Hz (dB)	60	52
		HFA-FOG (dB)	61	53
Reference test gain (dB)			53	42
Frequency range (Hz)			100-8900	100-7500
Telecoil output		1 mA/m field (1600 Hz) (dB SPL)	91	
		10 mA/m field (1600 Hz) (dB SPL)	111	
		HFA-SPLITS L/R (dB SPL)		100/100
Total harmonic distortion (Input 70 dB SPL)		500 Hz (%)	<9	<2
		800 Hz (%)	<6	<2
		1600 Hz (%)	<3	<2
Equivalent input noise level		Omni (dB SPL)	17	16
		Dir (dB SPL)	26	28
Battery			Lithium-ion	Lithium-ion
Expected operating time, hours ²				24


1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+AMD1:1994 but without influence of feedback.

2) Expected operating time for rechargeable battery depends on use pattern, active feature set, hearing loss, sound environment, battery age and use of wireless accessories.

		Ear Simulator	2CC Coupler	
		Measured according to IEC 60118-0:1983/AMD1:1994, IEC 60118-0:2015, IEC 60118-1:1995+AMD1:1998 CSV and IEC 60318-4:2010	Measured according to ANSI S3.22-2014, IEC 60118-0:2015 and IEC 60318-5:2006	
 100 Power receiver mold, Bass & Power dome Technical information Omnidirectional mode is used unless otherwise stated. Warning to the hearing aid dispenser The maximum output capability of the hearing aid may exceed 132 dB SPL (IEC 711). Special care should be exercised in selecting and fitting the hearing aid, as there may be risk of impairing the remaining hearing of the hearing aid user. — Acoustic input: 60 dB SPL - - - Magnetic input: 31.6 mA/m		OSPL90 	OSPL90  Full-on Gain  Frequency Response 	
OSPL90		Peak (dB SPL)	132	124
		1600 Hz (dB SPL)	130	122
		HFA-OSPL90 (dB SPL)	127	120
Full-on Gain ¹		Peak (dB)	66	57
		1600 Hz (dB)	60	52
		HFA-FOG (dB)	61	53
Reference test gain (dB)			53	42
Frequency range (Hz)			100-7500	100-7500
Telecoil output		1 mA/m field (1600 Hz) (dB SPL)	91	
		10 mA/m field (1600 Hz) (dB SPL)	111	
		HFA-SPLITS L/R (dB SPL)		100/100
Total harmonic distortion (Input 70 dB SPL)		500 Hz (%)	<9	<2
		800 Hz (%)	<6	<2
		1600 Hz (%)	<3	<2
Equivalent input noise level		Omni (dB SPL)	17	17
		Dir (dB SPL)	26	29
Battery			Lithium-ion	Lithium-ion
Expected operating time, hours ²			24	








1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+AMD1:1994 but without influence of feedback.

2) Expected operating time for rechargeable battery depends on use pattern, active feature set, hearing loss, sound environment, battery age and use of wireless accessories.



1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+AMD1:1994 but without influence of feedback.

2) Expected operating time for rechargeable battery depends on use pattern, active feature set, hearing loss, sound environment, battery age and use of wireless accessories.

		Ear Simulator Measured according to IEC 60118-0:1983/AMD1:1994, IEC 60118-0:2015, IEC 60118-1:1995+AMD1:1998 CSV and IEC 60318-4:2010	2CC Coupler Measured according to ANSI S3.22-2014, IEC 60118-0:2015 and IEC 60318-5:2006
 105 Power receiver mold		OSPL90   Full-on Gain   Frequency Response   — Acoustic input: 60 dB SPL - - - Magnetic input: 31.6 mA/m	
	Peak (dB SPL)	135	127
OSPL90	1600 Hz (dB SPL)	133	126
	HFA-OSPL90 (dB SPL)	131	123
	Peak (dB)	72	64
Full-on Gain ¹	1600 Hz (dB)	66	59
	HFA-FOG (dB)	65	58
	Reference test gain (dB)	58	47
	Frequency range (Hz)	100-7500	100-7500
	1 mA/m field (1600 Hz) (dB SPL)	96	
Telecoil output	10 mA/m field (1600 Hz) (dB SPL)	116	
	HFA-SPLITS L/R (dB SPL)		104/104
	500 Hz (%)	<2	<2
Total harmonic distortion (Input 70 dB SPL)	800 Hz (%)	<2	<2
	1600 Hz (%)	<4	<2
	Omni (dB SPL)	16	16
Equivalent input noise level	Dir (dB SPL)	25	28
Battery		Lithium-ion	Lithium-ion
Expected operating time, hours ²		24	

1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+AMD1:1994 but without influence of feedback.

2) Expected operating time for rechargeable battery depends on use pattern, active feature set, hearing loss, sound environment, battery age and use of wireless accessories.

[illegible]

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Headquarters
Oticon A/S
Kongebakken 9
DK-2765 Smørum
Denmark



SBØ Hearing A/S
Kongebakken 9
DK-2765 Smørum
Denmark

26341205 / 2022.11.07 / v1