



OtoNova

USE	APPLICATION PACKS	SCREENER	CLINICAL
	NEWBORN SCREENING PAEDIATRIC ASSESSMENT CLINICAL AUDIOLOGY RESEARCH		
DP OAE		DPOAE only	DPOAE
	TEST CONFIGURABILITY	9 Parameters, 1 Protocol	27 Parameters, Unlimited customisable Protocols
	PRESET FREQUENCY (kHz)	Up to 8kHz, 7 frequencies	Up to 8kHz, 16 frequencies
	FREQUENCY RANGE (kHz)	1 - 8kHz	0.75 - 8kHz
	FREQUENCY RESOLUTION	Half octave	Up to 100 pts/octave with DP Zoom
	FREQUENCY RATIO	1.22 (Fixed)	1.22 (Fixed)
	STIMULUS LEVELS (dBSPL)	60/55 or 70/70 (6 Presets)	60/50 to 70/70 (6 Presets) Custom 40 to 75 (Δ1dB)
	DISPLAY	Graphic progress indicators / Numeric / Paediatric distraction / Line Graph / SNR Histogram	Graphic progress indicators / Numeric / Paediatric distraction / Line Graph / SNR Histogram
	PASS CRITERIA	2, 3, 4 or 5 frequencies	Up to 16 freq. or 10-100% freq with DP Zoom
	PASS SNR	6dB (Fixed), 2 noise criteria levels	6dB, configurable up to 14dB (Δ1dB)
	FEATURES	Rapid DPOAE screening mode Paediatric distraction display Easy-to-use pass/refer result Automatic test progression and interpretation Large searchable database	Extensive configurability incl.12 frequency DPgram Hi-res DP-Zoom to examine DPgram features Programmable clinical assessment decisions Includes all SCREENER functions Easy switching between modes
TE OAE		TEOAE only ADD TE to DPOAE	INCLUDED
	TEST CONFIGURABILITY	8 Parameters, 1 Protocol	32 Parameters, Unlimited customisable Protocols
	FREQUENCY RANGE (kHz)	1, 1.5, 2, 3, 4 kHz	1, 1.5, 2, 3, 4, 6 kHz
	FREQUENCY RESOLUTION	Half octave bandwidth	Half octave bandwidth
	STIMULUS LEVELS (dBSPL _{pe})	Fixed 84 dBSPL _{pe}	40 to 90 (Δ1dB)
	DISPLAY	Graphic progress indicators / Numeric / Paediatric distraction / Histogram / Waveform	Graphic progress indicators / Numeric / Paediatric distraction / Histogram / Waveform / Spectrum
	PASS CRITERIA	2, 3, or 4 half octave bands or narrow bands	1 to 6 half octave bands, settable group SNR
	PASS SNR (dB)	6dB (Fixed)	Configurable up to 14dB (Δ1dB) for automation
	FEATURES	Gold standard TEOAE Quickscreen, ideal for infant screening Easy-to-use fast pass/refer result Automatic test progress and interpretation Large searchable database	1-6kHz range, half octave resolution Hi-res TEOAE spectrum Programmable clinical assessment decisions. Custom Norms display Includes all SCREENER functions Easy switching between modes
	AABR		AABR only ADD AABR to OAE
TEST CONFIGURABILITY		7 Parameters, 1 Protocol	14 Parameters, Unlimited Protocols
STIMULUS LEVELS (dBnHL)		30 to 45 (Δ1dB) 40 with earcup	20 to 45 (Δ1dB)* 40 with earcup
STIMULUS TYPE		Chirp	Chirp and Click
FEATURES		Rapid AABR screening mode Easy-to-use automatic pass/refer result Chirp stimulation with weighted averaging Superior pretest environment checks Powerline and individual impedance analysis Complete electrical isolation Earcup or insert stimulator	Chirp or click stimulation with weighted averaging **60 dBnHL~ for training insert only Includes easy access to OtoNova AABR SCREENER mode features
ADDITIONAL OPTIONS			
Wireless Module		-	-
Docking Station		-	-
Integral Barcode / QR Scanner		-	-
Add Binaural Testing		OPTION	OPTION
Wireless Printer	-	-	
Add Carry Case from	OPTION	OPTION	
* subject to qty			

NEW



HANDHELD

Otoport

Features

Common to all Otoport models

- Light weight, easy-to-use screening or clinical device
- Re-configurable for DPOAE and/or TEAOE and/or AABR
- Large database with extensive patient data fields
- Programmable work-list
- Large battery capacity - over 250 tests per battery charge
- Create multiple user profiles with Admin management and security
- On-screen notifications probe fit, noise levels and stimulus quality
- Clear test progression graphics and notification
- Each test is saved with the actual psychophysical response data
- Transfer data to the Otolinik pc App wirelessly or via USB
- Automated self-test system when powering on Otoport
- Probe function test and ID tracking
- Customisable QA calibration test area & history view
- Optional Built in QR and barcode scanner
- Print wirelessly to the Otodynamics optional mini printer

Otoport



Otoport

LITE	FLEXI	ADVANCE
DPOAE only	DPOAE	DPOAE
11 Parameters, 1 Protocol 2, 3, 4, 6 kHz, 4 frequencies 2 - 6 kHz Half octave 1.22 (Fixed) 60/55 or 65/50 or 65/55 Graphic progress indicators / Numeric / Paediatric distraction / Line Graph Pass on 2, 3 or 4 frequencies 6dB (Fixed) Rapid DPOAE screening mode Interactive Paediatric distraction display Easy-to-use pass/refer result Automatic test progression and interpretation	20 Parameters, 4 Protocols Up to 8 kHz, 13 frequencies 1 - 8 kHz Up to 4 pts/octave 1.22 (Fixed) 60/50 to 70/70 (6 Presets) Graphic progress indicators / Numeric / Paediatric distraction / Line Graph / Histogram 2 to 6 frequencies, 12 frequency DP 6dB configurable up to 14dB (Δ 1dB) Extensive configurability including 12 frequency DPOAE Configurable norms Includes all LITE functions Easy switching between applications	24 Parameters, 8 Protocols Up to 8kHz, 16 frequencies (audiometric or extract values) 0.75 to 8 kHz Up to 100 pts/octave with DP Zoom 1.22 (Fixed) 60/50 to 70/70 (6 Presets) Custom 40 to 75 (Δ 1dB) Graphic progress indicators / Numeric / Paediatric distraction / Line Graph / Histogram 1 to 6 DP or 12 freq DP or 10-100% freq. with DP zoom 6dB configurable up to 14dB (Δ 1dB) DP Zoom to 100pts/octave DP Growth I/O with scissor mode Access to selected research function for IRB approved projects. Multi test macro (4x protocols) Includes all FLEXI functions
TEOAE only ADD TEOAE to DP	TEOAE INCLUDED	TEOAE INCLUDED
14 Parameters, 1 Protocol 1, 1.5, 2, 3, 4 kHz Half octave bandwidth Fixed 84 dB SPLpe Graphic progress indicators / Numeric / Paediatric distraction / Histogram Customisable SNR in any of 2, 3, or 4 half octave bands, or wide / narrow bands 6dB (Fixed) Gold standard TEOAE Quickscreen Interactive Paediatric distraction display Easy-to-use pass/refer result Automatic test progress and interpretaon	25 Parameters, 2 Protocols 1, 1.5, 2, 3, 4, 6 kHz Half octave bandwidth 60 to 90 dB SPLpe (Δ 1dB) Graphic progress indicators / Numeric / Paediatric Distraction / Histogram / TE Waveform Includes all LITE functions. Easy switching between applications with additional protocol.	28 Parameters, 2 Protocols 1, 1.5, 2, 3, 4, 6 kHz Half octave bandwidth 60 to 90 dB SPLpe (Δ 1dB) Graphic progress indicators / Numeric / Paediatric Distraction / Histogram / TE Waveform / high res TE spectrum / Stimulus Waveform & Spectrum Customisable SNR in any of 1 to 5 half octave bands 6dB recommended variable up to 14dB (Δ 1dB) Access to Otoptop FLEXI features
AABR only ADD AABR to OAE	ADD AABR to OAE	ADD AABR to OAE
8 Parameters, 1 Protocol 30 to 45 (Δ 1dB) (up to 40 with earcups) Chirp Superior pretest environment, with power-line & individual impedance analysis. Fast pass/refer result Choice in stim delivery ear cup/insert	11 Parameters, 2 Protocols 30 to 45 (Δ 1dB) (up to 40 with earcups)** Chirp Real time ABR waveform view ** 45 - 60 dBnHL insert only for training Includes all LITE functions. Easy switching between applications with additional protocol.	16 Parameters, 2 Protocols 30 to 45 (Δ 1dB) (up to 40 with earcups)** Chirp and Click Real time ABR waveform view **45 to 60 dBnHL insert only for training Extendable FSP for AABR waveform capture Easy access to all Otoptop FLEXI features
ADDITIONAL OPTIONS		
INCLUDED	INCLUDED	INCLUDED
INCLUDED	INCLUDED	INCLUDED
OPTION	OPTION	OPTION
OPTION	OPTION	OPTION
OPTION	OPTION	OPTION
OPTION	OPTION	OPTION

USE

DP OAF

ECAF

Features

Common to both
Otocheck LE models

- Easy-to-use pass/refer result
- 1 key press and go operation
- Minimal training required
- SNR & Signal to noise histogram display
- Interactive Paediatric distraction display
- Individual frequency progress display
- 250 tests per battery charge

ADDITIONAL OPTIONS	
Wireless Module	INCLUDED
Charging Cradle	INCLUDED
Wired Printer	OPTION
Add Carry Case from	OPTION

Otoport



Echoport



- Laptop/Tablet not included

Otodynamics Ltd.

30-38 Beaconsfield Road Hatfield Herts AL10 8BB UK
T. +44 (0) 1707 267540 E. info@otodynamics.com
www.otodynamics.com

(Product specification may differ)

Ref: F-008-19ROW REV005